

# Air Townsend & Air Lite

## Ligament Bracing Solutions

### Routinely Prescribed For...

Patients who desire or require maximum control after a ligament injury and/or reconstructive knee surgery.

### This Product Can Be Ordered...

- Custom from a cast, digital scan or KneePix

### Offered With These Standard Features...

- Custom Molded Laminated Carbon Graphite Shells
- Patented Townsend Motion TM5+ Hinges
- Extension Stop Kit
- Patented Synergistic Suspension Strap
- Anti-Migration Strap Padding
- Anti-Rotation Tibia Shell Bolster (Air Townsend)

### Available Options Include...


- ACL Or Combined Instabilities
- Virtually Any Shell Length Over 13 Inches
- Hinge Material -- Aluminum, Titanium Or Stainless
- Extension Assist
- Special Trim Lines For Skiers, Rodeo, Hockey Goalies, And Full Figure (Air Townsend) Model For Soft Tissue Containment
- Flexion Stop Kit
- A Variety Of Colors Or Fabric Inlay



The Air Townsend is a traditional full shell design that provides maximum control and protection for ligament instabilities.



The Air Lite features an adjustable anterior tibia strap and narrower bands to reduce the weight and profile of the brace.



## Optional Shell Configuration

Optional single upright KAFO with leg extension and heel cup for rotary control or to eliminate any potential for brace migration.

# General Information

Patient's Last Name: \_\_\_\_\_

Patient's First Name: \_\_\_\_\_

Male  Female Age \_\_\_\_\_ Height \_\_\_\_\_ Weight \_\_\_\_\_

Work/Activities:  Activities of Daily Living

Non-Contact Sports  Contact Sports

Leg:  Left  Right

Ligament:  ACL  PCL  LCL  MCL

Meniscus Damage:  Medial  Lateral

Surgeries (type/date): \_\_\_\_\_

Prescribing Physician: \_\_\_\_\_

If there is a question about this order, who should we contact?

Name: \_\_\_\_\_ Phone: \_\_\_\_\_

Email: \_\_\_\_\_

BILLING: P.O. Number \_\_\_\_\_

Townsend Account Number: \_\_\_\_\_

Bill To: \_\_\_\_\_

Address: \_\_\_\_\_

City: \_\_\_\_\_

State: \_\_\_\_\_ Zip Code: \_\_\_\_\_ Country: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Ship To: \_\_\_\_\_

Address: \_\_\_\_\_

City: \_\_\_\_\_

State: \_\_\_\_\_ Zip Code: \_\_\_\_\_ Country: \_\_\_\_\_

Phone: \_\_\_\_\_ Attention: \_\_\_\_\_

**Fit Date: If known, please indicate the date you are scheduled to fit the patient:** \_\_\_\_\_

## Shipping Preference:

Ground  2-Day P.M.  2-Day A.M.

Next Day P.M.  Next Day A.M.

(If no preference is indicated, this order will be shipped 2 Day P.M.)

**Note:** We do not ship new or repaired products directly to patients.

# Air Townsend & Air Lite Ligament Knee Braces

Casts: 18-20 inch length; full extension; non-weight bearing, foot dorsiflexed; quad relaxed; landmarks indicated; cut off back of the leg.

**\*Indicates additional charges apply**

Air Townsend  Air Lite

## Model

ACL  Combined Instabilities  
 Thigh Strap\*  Graphite Band\*

## Thigh Shell Length

7 Inch  8 Inch  Other \_\_\_\_\_"

## Thigh Band Width

1.5 Inches  2 Inches  2.5 Inches

## Tibia Shell Length

7 Inch  8 Inch  9 Inch

Single Strut KAFO With Heel Cup

(Must Complete Additional Form For Brace Extension)

## Hinges

TM5+ Hinges -- Includes extension stop kit

Optional flexion stop kit\*

Add optional extension assist bands/posts\*

## Hinge Material

6061 Aluminum (standard, if no hinge material is indicated)

Stainless Steel\*  Titanium\*

Original Hinges\* -- Titanium (standard) includes extension stop kit

Condylar Pads:  No  Bi-Lateral  Medial  Lateral

Optional flexion stop kit\*

Add optional extension assist bands/posts\*

**Tibia Shell Anti-Rotation Bolster** (Air Townsend only)

Standard  Minimal  Soft  No Bolster

**Special Trim Lines** (Air Townsend only)

Full Figure: For additional soft tissue containment

No tibia shell hole (Townsend Original)

Ski boot cut: 1/2 inch notched in distal tibia shell

Rodeo: Medial thigh cut-out for rider comfort

Customized Shell Design (include instructions/drawing)

**Synergistic Suspension Strap Attachment**

Standard attachment (Lateral end recessed inside tibia shell)

Lateral end attached to outer shell (for prominent fibular head)

## Color/Fabric Inlay

Black  Beige  Gray  Red  Navy Blue

Royal Blue  Green  Burgundy

Clear Graphite  Sheer Red\*

Sheer Teal\*  Sheer Purple\*

Fabric -1 yard from patient\*  U.S.A. Flag Fabric\*

Spooner Patella Stabilizing Attachment\*

## Brace Cover\*

Posterior Closure:  Black  Blue

Pull On

## Undersleeves\*

18" Cotton  18" Neoprene  22" Neoprene

## Thigh Sleeves\*

1/8 Atrophy Thigh Sleeve  1/16 Comfort Thigh Sleeve

## Cast Tape (select type, quality and tape width)

FAB SERIES (hybrid polyester) or  FG SERIES (fiber glass)

Quantity:  2 Rolls  4 Rolls  6 Rolls

Width:  2"  3"  4"  5"

**Special Instructions:** \_\_\_\_\_

## TOWNSEND'S SHIPPING DEPARTMENT USE ONLY

NEW BRACE  SERVICE -- Original Brace Returned?  Yes  No

RECEIVED \_\_\_\_\_ SHIPPED VIA \_\_\_\_\_

## Townsend Design

4615 Shepard St., Bakersfield, CA, 93313  
Phone: 800.432.3466 or 661.837.1795; Fax: 800.798.2722