

Walking Boots

Post Operative Bracing Solutions

Routinely Prescribed For...

Patients who require protection and/or stabilization to ambulate while they recover from foot, toe or ankle injuries and/or surgery. *This product is contraindicated for unstable fractures or fractures of the tibia or fibula.*

This Product Can Be Ordered...

- Standard sizes (based on patient's shoe size)
- TD plastic walkers are available S-XL
- TSB metal boots are sizes S,M,L

Offered With These Standard Features...

- Foot Shell
- Aluminum Uprights (that you can contour)
- Molded/Padded Insole
- Shock Absorbing Outer Sole Rocker
- Foam Bootie
- Non-Elastic Straps

Available Options Include...

- Plastic Or Metal Foot Shell
- Full Height Or Mid Calf Uprights
- ROM Hinges
- Standard Or Pneumatic (Air Bladder) Bootie
- Heel Wedges
- Indoor Sanitary Outer Sole Cover



TD Walker



TD Mid Calf Walker



TD ROM Walker

Townsend's Stability Boot Series provides a range of ordering options, including rigid metal shell walkers and economical plastic shell walking boots.



Townsend's top-of-the-line metal shell TSB features an adjustable, semi-rigid posterior shell. This feature allows you to align the midline of the leg with the center of the uprights to ensure the proper weight bearing line is achieved.



AIR Walkers feature a pneumatic bootie.



Accessories include a heel wedge kit and hygienic indoor sole cover.

Ordering Options

Townsend's plastic shell walkers are the lightest weight boots available and are a low profile, low cost solution.

Plastic shell boots are available with ROM hinges.

The uprights on Mid Calf Walkers are approximately 4 inches shorter than the uprights on full height walkers.

Pneumatic booties are an additional charge item that can be substituted for the standard bootie on all walker models.



www.townsenddesign.com

Toll Free: 800.432.3466 Ph: 661.837.1795 Fax: 661.837.0613

General Information

Ordered by : _____

Phone: _____

Patient's Name: _____

BILLING: P.O. Number _____

Townsend Account Number: _____

Bill To: _____

Address: _____

City: _____

State: _____ Zip Code: _____ Country: _____

Phone: _____ Fax: _____

Ship To: _____

Address: _____

City: _____

State: _____ Zip Code: _____ Country: _____

Phone: _____ Attention: _____

Shipping Preference:

Ground 2-Day P.M. 2-Day A.M.

Next Day P.M. Next Day A.M.

(If no preference is indicated, this order will be shipped 2 Day P.M.)

Note: We do not ship new or repaired products directly to patients.

Sizing Reference

Small: Men's Shoe Size 3.5-7; Women's Shoe Size 4.5-8

Medium: Men's Shoe Size 7-10; Women's Shoe Size 8-11

Large: Men's Shoe Size 10-13; Women's Shoe Size 11-15

Extra Large: Men's Shoe Size 13-16; Women's Shoe Size 15+

Note: Upright heights vary based on the boot model and size. The uprights on Mid-Calf walkers are approximately four inches shorter than the uprights on full height walkers

Ordering Instructions

Please complete and fax this form to 800.798.2722 (24-hours a day). If you are calling in your order, this form indicates the options and information that will be required by our staff. For phone orders, please call 800.700.2722 between 7:00 a.m. and 4:00 p.m. (PST).

Townsend Walking Boots

Plastic Shell Walkers

TD Walker -- Quantity: S___ M___ L___ XL___
with plastic shells and aluminum struts

TD MC Walker -- Quantity: S___ M___ L___ XL___
with shorter mid-calf struts

TD ROM Walker -- Quantity: S___ M___ L___ XL___
with adjustable ROM hinges

TD AIR Walker -- Quantity: S___ M___ L___ XL___
with pneumatic bootie

TD MC AIR Walker -- Quantity: S___ M___ L___ XL___
with mid-calf struts and pneumatic bootie

TD AIR ROM Walker -- Quantity: S___ M___ L___ XL___
with ROM hinges and pneumatic bootie

Metal Shell Walkers

TSB -- Quantity: S___ M___ L___
with metal shells and StableLock heel positioning shell

TSB-MC -- Quantity: S___ M___ L___
with shorter struts and StableLock heel positioning shell

TSB Air -- Quantity: S___ M___ L___
with metal shells and StableLock heel positioning shell

TSB-MC Air -- Quantity: S___ M___ L___
metal shells with shorter struts and StableLock heel positioning shell

Please Indicate Accessories

Walker Wedges* -- Quantity: ___ S___ M___ L___ XL___
angled heel wedge insert kit

Indoor Walker Sole Cover* -- Quantity: ___
keeps floors/carpets clean

Replacement Liners

Standard Mid-Calf Liners -- Quantity: S___ M___ L___ XL___

Standard Tall Liners -- Quantity: S___ M___ L___ XL___

Pneumatic Mid-Calf Liners -- Quantity: S___ M___ L___ XL___

Pneumatic Tall Liners -- Quantity: S___ M___ L___ XL___



TD MC Walker

Plastic foot shell and aluminum uprights



TD Walker

Plastic foot shell and aluminum uprights



TD ROM Walker

Plastic foot shell, ROM hinges, and aluminum uprights



TSB

Metal foot shell and aluminum uprights

Townsend Design

4615 Shepard St., Bakersfield, CA, 93313

Phone: 800.432.3466 or 661.837.1795; Fax: 800.798.2722