

Reliever Air & Reliever Air Lite

OA Bracing Solutions

Routinely Prescribed For...

Patients who have moderate to severe medial or lateral unicompartement osteoarthritis and/or combination OA and ligament instability. Ideal for patients requiring soft tissue containment and maximum unloading.

This Product Can Be Ordered...

- Custom from a cast, digital scan or KneePix

Offered With These Standard Features...

- Custom Molded Laminated Carbon Graphite Shells
- Patented Townsend Motion TM5+ Hinges
- Extension Stop Kit
- Patented Loadshifter
- Patented Synergistic Suspension Strap
- Anti-Migration Strap Padding
- Anti-Rotation Tibia Shell Bolster (Reliever Air)

Available Options Include...

- Virtually Any Shell Length Over 13 Inches
- Hinge Material -- Aluminum, Titanium Or Stainless
- Extension Assist
- Full Figure Model For Soft Tissue Containment
- Flexion Stop Kit
- A Variety Of Colors Or Fabric Inlay



The Reliever Air is a traditional full shell design that provides maximum unloading, soft tissue containment and support for ligament instabilities.



The Lite model features an adjustable anterior tibia strap and modified trim lines that reduce the weight of the brace.

Corrective Force Technology



TOWNSEND DESIGN
LOADSHIFTER
RELIEF MECHANISM

The Reliever Air and Reliever Air Lite feature Townsend's patent Loadshifter for adjusting the angle of the thigh shell to increase the base level of corrective force built into the brace.



Shell Options

Optional single upright KAFO with leg extension and heel cup for rotary control or suspension.

 **Townsend**

www.townsenddesign.com

Toll Free: 800.432.3466 Ph: 661.837.1795 Fax: 661.837.0613

General Information

Patient's Last Name: _____

Patient's First Name: _____

Male Female Age ____ Height ____ Weight ____

Work/Activities: Activities of Daily Living

Non-Contact Sports Contact Sports

Leg: Left Right

Ligament: ACL PCL LCL MCL

Meniscus Damage: Medial Lateral

Surgeries (type/date): _____

Prescribing Physician: _____

If there is a question about this order, who should we contact?

Name: _____ Phone: _____

Email: _____

BILLING: P.O. Number _____

Townsend Account Number: _____

Bill To: _____

Address: _____

City: _____

State: _____ Zip Code: _____ Country: _____

Phone: _____ Fax: _____

Ship To: _____

Address: _____

City: _____

State: _____ Zip Code: _____ Country: _____

Phone: _____ Attention: _____

Fit Date: If known, please indicate the date you are scheduled to fit the patient: _____

Shipping Preference:

Ground 2-Day P.M. 2-Day A.M.

Next Day P.M. Next Day A.M.

(If no preference is indicated, this order will be shipped 2 Day P.M.)

Note: We do not ship new or repaired products directly to patients.

Full Shell Reliever Series

OA Unloading Knee Braces

Casts: 18-20 inch length; full extension; non-weight bearing, foot dorsiflexed; quad relaxed; landmarks indicated; cut off back of the leg.

*Indicates additional charges apply

Compartment Unload Medial Compartment
 Unload Lateral Compartment

Check Box Below To Select Model



Thigh Shell

7 Inch 8 Inch 9 Inch

Tibia Shell

7 Inch 8 Inch 9 Inch Other: _____

Single Strut KAFO With Heel Cup

(Must Complete Additional Form For Brace Extension)

Brace is fabricated WITH LOADSHIFTER (no charge)

The LOADSHIFTER can be used to increase correction by shifting the thigh shell angle. ONLY offered with TM5+ Aluminum Hinges.

TM5+ Hinges -- Includes extension stop kit

Optional flexion stop kit*

Hinge Material

6061 Aluminum (required for LOADSHIFTER)

Stainless Steel* (brace will be fabricated without LOADSHIFTER)

Titanium* (brace will be fabricated without LOADSHIFTER)

Add optional extension assist bands/posts*

Tibia Shell Anti-Rotation Bolster (Reliever Air & Reliever only)

Standard Minimal Soft No Bolster

Special Trim Lines

Full Figure (Reliever Air Only): Helps contain soft tissue

Customized Shell Design (include instructions/drawing)

Synergistic Suspension Strap Attachment

Standard attachment (Lateral end recessed inside tibia shell)

Lateral end attached to outer shell (for prominent fibular head)

Color/Fabric Inlay

Black Beige Gray Red Navy Blue
 Royal Blue Green Burgundy

Clear Graphite Sheer Red*

Sheer Teal* Sheer Purple*

Fabric -1 yard from patient* U.S.A. Flag Fabric*

Spooner Patella Stabilizing Attachment*

Brace Cover*

Posterior Closure: Black Blue

Pull On

Undersleeves*

18" Cotton 18" Neoprene 22" Neoprene

Thigh Sleeves*

1/8 Atrophy Thigh Sleeve 1/16 Comfort Thigh Sleeve

Cast Tape (select type, quality and tape width)

FAB SERIES (hybrid polyester) or FG SERIES (fiber glass)

Quantity: 2 Rolls 4 Rolls 6 Rolls

Width: 2" 3" 4" 5"

Special Instructions: _____

TOWNSEND'S SHIPPING DEPARTMENT USE ONLY

NEW BRACE SERVICE -- Original Brace Returned? Yes No

RECEIVED _____ SHIPPED VIA _____

Townsend Design

4615 Shepard St., Bakersfield, CA, 93313
Phone: 800.432.3466 or 661.837.1795; Fax: 800.798.2722