

# Premier Series

## Ligament Bracing Solutions

### Routinely Prescribed For...

Patients who have experienced ligament injuries and/or reconstructive knee surgery.

### This Product Can Be Ordered...

- Custom from a cast, digital scan, tracing, CAS or KneePix

### Offered With These Standard Features...

- Custom Molded Solid Core Carbon Graphite Shells
- Patented Townsend Motion TM5+ Hinges
- Extension Stop Kit
- Patented Synergistic Suspension Strap
- Anti-Migration Strap Padding
- Rotation Control Tibial Shell Contour
- High Gloss Paint Or Powdercoat Finish

### Available Options Include...

- ACL Or Combined Instabilities
- 13, 14, 15 And 16 Inch Shell Length
- Hinge Material -- Aluminum, Titanium Or Stainless
- Compression/Suspension Package
- Flexion Stop Kit
- A Variety Of Colors



Townsend's most popular custom ligament knee brace, the Premier features a lightweight, low profile design that achieves excellent control and patient compliance. While anterior shells are recommended, an optional posterior tibia band is available.

## Optional Shell Configurations

The Premier is available in a variety of custom shell configurations.



Optional Maximum Control (MC) thigh has increased surface contact for added stabilization.



Optional double band tibia provides additional support and functional control.



Optional posterior tibia shell.



Optional single upright KAFO with leg extension and heel cup for rotary control or to eliminate any potential for brace migration.



[www.townsenddesign.com](http://www.townsenddesign.com)

Toll Free: 800.432.3466 Ph: 661.837.1795 Fax: 661.837.0613

# General Information

Patient's Last Name: \_\_\_\_\_

Patient's First Name: \_\_\_\_\_

Male  Female Age \_\_\_\_ Height \_\_\_\_ Weight \_\_\_\_

Activities:  Activities of Daily Living

Non-Contact Sports  Contact Sports

Leg:  Left  Right

Ligament:  ACL  PCL  LCL  MCL

Meniscus Damage:  Medial  Lateral

Surgeries (type/date): \_\_\_\_\_

Prescribing Physician: \_\_\_\_\_

If there is a question about this order, who should we contact?

Name: \_\_\_\_\_ Phone: \_\_\_\_\_

Email: \_\_\_\_\_

BILLING: P.O. Number \_\_\_\_\_

Townsend Account Number: \_\_\_\_\_

Bill To: \_\_\_\_\_

Address: \_\_\_\_\_

City: \_\_\_\_\_

State: \_\_\_\_\_ Zip Code: \_\_\_\_\_ Country: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Ship To: \_\_\_\_\_

Address: \_\_\_\_\_

City: \_\_\_\_\_

State: \_\_\_\_\_ Zip Code: \_\_\_\_\_ Country: \_\_\_\_\_

Phone: \_\_\_\_\_ Attention: \_\_\_\_\_

**Fit Date: If known, please indicate the date you are scheduled to fit the patient:** \_\_\_\_\_

### Shipping Preference:

Ground  2-Day P.M.  2-Day A.M.

Next Day P.M.  Next Day A.M.

(If no preference is indicated, this order will be shipped 2 Day P.M.)

**Note:** We do not ship new or repaired products directly to patients.

# Premier Ligament Knee Braces

Casts: 18-20 inch length; full extension; non-weight bearing, foot dorsiflexed; quad relaxed; landmarks indicated; cut off back of the leg.

**\*Indicates additional charges apply**


Model  ACL  
 Combined Instabilities (choose strap or band)  
 PCL Strap\*  PCL Rigid Band\*

Thigh Shell Length  7 Inch  8 Inch

Tibia Shell Length  6 Inch  7 Inch  8 Inch

### Thigh

A: Standard  
 B: Maximum Control (MC)



### Tibia

C: Anterior Single Band  
 D: Posterior Single Band  
 E: Double Band\* (7" or 8")  
 Single Strut KAFO (Must Complete Additional Form For Brace Extension)



### TM5+ Hinges -- Includes extension stop kit

Optional Flexion Stop Kit\*  
 Add optional extension assist bands/posts\*

### TM5+ Hinge Material

6061 Aluminum  Stainless Steel\*  Titanium\*

### Finish and Color

Powdercoat Finish  
 Black  Antique Pewter  
 Royal Blue  Burgundy

### High Gloss Paint Finish

Black  Royal Blue  Burnt Orange  
 Dark Violet  Emerald Green  Steel Blue  
 Quicksilver  Indy Yellow  Burgundy  
 White  Beige

Custom Paint Finish\* -- Indicate Custom Paint # \_\_\_\_\_

### Options

C/S Package\*  
Compliance/Suspension semi rigid padded inserts added to thigh shell to increase dynamic compression and enhance suspension.  
 No wraparound attachment of Synergistic Suspension Strap (recommended if patient has a prominent fibular head)

Spooner Patella Stabilizing Attachment\*

### Brace Cover\*

Posterior Closure:  Black  Blue  
 Pull On

### Undersleeves\*

18" Cotton  18" Neoprene  22" Neoprene

### Thigh Sleeves\*

1/8 Atrophy Thigh Sleeve  1/16 Comfort Thigh Sleeve

### Cast Tape (select type, quality and tape width)

FAB SERIES (hybrid polyester) or  FG SERIES (fiber glass)

Quantity:  2 Rolls  4 Rolls  6 Rolls

Width:  2"  3"  4"  5"

### TOWNSEND'S SHIPPING DEPARTMENT USE ONLY

NEW BRACE  SERVICE -- Original Brace Returned?  Yes  No

RECEIVED \_\_\_\_\_ SHIPPED VIA \_\_\_\_\_

## Townsend Design

4615 Shepard St., Bakersfield, CA, 93313  
Phone: 800.432.3466 or 661.837.1795; Fax: 800.798.2722

Special Instructions: \_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_