

# Premier Series

## OA Bracing Solutions

### Routinely Prescribed For...

Patients who have moderate to severe medial or lateral unicompartement osteoarthritis. The double upright model is recommended for combination OA and ligament instability. The single hinge PR1 is ideal for lateral OA.

### This Product Can Be Ordered...

### Offered With These Standard Features...

- Custom Molded Solid Core Carbon Graphite Shells
- Patented Townsend Motion TM5+ Hinges (double upright)
- Patented Townsend Motion TM6 Hinge (single upright)
- Extension Stop Kit
- Patented Synergistic Suspension Strap
- Anti-Migration Strap Padding
- Rotation Control Tibial Shell Contour
- High Gloss Paint Or Powdercoat Finish

### Available Options Include...

- 14, 15 And 16 Inch Shell Length
- Hinge Material -- Aluminum, Titanium Or Stainless
- Extension Assist (double upright model)
- Compression/Suspension Package (double upright model)
- Flexion Stop Kit
- A Variety Of Colors



Premier Reliever features a lightweight, low profile design that provides comfortable unloading and pain relief. This picture shows the optional C/S Package that enhances suspension.




The new PR1 is a single upright model that is ideal for lateral OA (no medial hinge) and patients who prefer an ultra light weight brace.

## Corrective Force Technology



The Premier Reliever features Townsend's patent **Loadshifter** for adjusting the angle of the thigh shell to increase the base correction built into the brace.



The PR1 features an adjustable force buttress that allows the patient to optimize counter force pressure. The buttress also improves brace suspension.

## Shell Options



Optional double band tibia for dual hinge Premier Reliever braces.



Optional single upright KAFO with leg extension and heel cup for rotary control or suspension.

 **Townsend**

[www.townsenddesign.com](http://www.townsenddesign.com)

Toll Free: 800.432.3466 Ph: 661.837.1795 Fax: 661.837.0613

# General Information

Patient's Last Name: \_\_\_\_\_

Patient's First Name: \_\_\_\_\_

Male  Female Age \_\_\_\_ Height \_\_\_\_ Weight \_\_\_\_

Activities:  Activities of Daily Living

Non-Contact Sports  Contact Sports

Leg:  Left  Right

Ligament:  ACL  PCL  LCL  MCL

Meniscus Damage:  Medial  Lateral

Surgeries (type/date): \_\_\_\_\_

Prescribing Physician: \_\_\_\_\_

If there is a question about this order, who should we contact?

Name: \_\_\_\_\_ Phone: \_\_\_\_\_

Email: \_\_\_\_\_

BILLING: P.O. Number \_\_\_\_\_

Townsend Account Number: \_\_\_\_\_

Bill To: \_\_\_\_\_

Address: \_\_\_\_\_

City: \_\_\_\_\_

State: \_\_\_\_\_ Zip Code: \_\_\_\_\_ Country: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Ship To: \_\_\_\_\_

Address: \_\_\_\_\_

City: \_\_\_\_\_

State: \_\_\_\_\_ Zip Code: \_\_\_\_\_ Country: \_\_\_\_\_

Phone: \_\_\_\_\_ Attention: \_\_\_\_\_

**Fit Date: If known, please indicate the date you are scheduled to fit the patient:** \_\_\_\_\_

### Shipping Preference:

Ground  2-Day P.M.  2-Day A.M.

Next Day P.M.  Next Day A.M.

(If no preference is indicated, this order will be shipped 2 Day P.M.)

**Note:** We do not ship new or repaired products directly to patients.

# Premier Reliever & Premier Reliever1

OA Unloading Knee Braces

Casts: 18-20 inch length; full extension; non-weight bearing, foot dorsiflexed; quad relaxed; landmarks indicated; cut off back of the leg.

▶ Indicates this option is not offered on Reliever1 model  
\*Indicates additional charges apply

### Model

Premier Reliever (dual hinge knee brace)

Premier Reliever1 (single hinge knee brace)

### Compartment

Medial Compartment  Lateral Compartment

### Thigh Shell Length

7 Inch  8 Inch

### Tibia Shell Length

7 Inch  8 Inch

### Tibia Shell

A: Single Band

▶  B: Anterior Double Band\*

C: Posterior Band

Single Strut KAFO With Heel Cup

(Must Complete Additional Form For Brace Extension)



**Premier Reliever** is fabricated with TM5+ Hinges and includes an adjustable extension stop kit

Optional Flexion Stop Kit\*

Add optional extension assist bands/posts\*

**Hinge Material Options** (dual hinge Premier Reliever only)

6061 Aluminum (required for adjustable correction LOADSHIFTER)

Stainless Steel (brace will be fabricated without LOADSHIFTER)

Titanium (brace will be fabricated without LOADSHIFTER)

**Premier Reliever1** is fabricated with an aluminum TM6 Hinge and includes an adjustable extension stop kit.

Optional Flexion Stop Kit

### Finish and Color

Powdercoat Finish

Black  Antique Pewter (Silver)

Royal Blue  Burgundy

High Gloss Paint Finish

Black  Royal Blue  Burnt Orange  Dark Violet

Emerald Green  Steel Blue  Quicksilver

Indy Yellow  Burgundy  White  Beige

Custom Paint Finish\* Indicate Custom Paint # \_\_\_\_\_

### Options

▶  **C/S Package\*** (only offered on dual hinge Premier Reliever)  
Compliance/Suspension semi rigid padded inserts added to thigh shell to increase dynamic compression and enhance suspension.

### Synergistic Suspension Strap

No wraparound attachment of Synergistic Suspension Strap (recommended if patient has a prominent fibular head)

### Accessories

Spooner Patella Stabilizing Attachment\*

### Brace Cover\*

Posterior Closure:  Black  Blue

Pull On

### Undersleeves\*

18" Cotton  18" Neoprene  22" Neoprene

### Thigh Sleeves\*

1/8 Atrophy Thigh Sleeve  1/16 Comfort Thigh Sleeve

### Cast Tape (select type, quality and tape width)

FAB SERIES (hybrid polyester) or  FG SERIES (fiber glass)

Quantity:  2 Rolls  4 Rolls  6 Rolls

Width:  2"  3"  4"  5"

### TOWNSEND'S SHIPPING DEPARTMENT USE ONLY

NEW BRACE  SERVICE -- Original Brace Returned?  Yes  No

RECEIVED \_\_\_\_\_ SHIPPED VIA \_\_\_\_\_

## Townsend Design

4615 Shepard St., Bakersfield, CA, 93313  
Phone: 800.432.3466 or 661.837.1795; Fax: 800.798.2722

Special Instructions: \_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_