




ROM-SS
(Breeze Pads)



The new ROM-SS (Suspension Solution) is the first post-op knee brace specifically designed to reduce brace migration. Featuring Townsend's C/S Package and semi-malleable aluminum cuffs that can be contoured to the sides of the leg, this brace will suspend better than any traditional post-operative knee brace. ROM-SS braces can be ordered with cooler Breeze pads or Wrap-Around pads. The ROM-SS also features a patient-friendly drop lock mechanism and tool-free length adjustment -- the standard model extends from 17-22 inches and the longer model extends from 19-24 inches. Rotating pull arms, under the cap, make it easy to set extension and flexion stops,



ROM6

Additional Models

ROM4: 4 strap, 17 inches, Full pads

Air ROM: 4 strap, 17 inches, Air pads

ROM6: 6 Strap, 23 inches, removable bar ends, Full pads

Air ROM-L: 4 strap, 23 inches, removable bar ends, Air pads

ROM6-L: 6 strap, 25 inches, removable bar ends, Full pads

ROM Knee Braces

Post Operative Bracing Solutions

Routinely Prescribed For...

Patients who require post-operative fixation and/or range of motion control for post-injury stabilization.

This Product Can Be Ordered...

- Off-The-Shelf in a variety of models & lengths

Offered With These Standard Features...

- Rigid Aluminum Bars
- Hinges With Adjustable Extension & Flexion Stops
- Drop Lock Mechanism

Available Options Include...

- Bar Length
- Shell Pads (Full, Air, Or Breeze)
- Straps With Velcro Or Buckle Closures

www.townsenddesign.com

Toll Free: 800.432.3466 Ph: 661.837.1795 Fax: 661.837.0613

 **Townsend**

General Information

Ordered by : _____

Phone: _____

Patient's Name: _____

BILLING: P.O. Number _____

Townsend Account Number: _____

Bill To: _____

Address: _____

City: _____

State: _____ Zip Code: _____ Country: _____

Phone: _____ Fax: _____

Ship To: _____

Address: _____

City: _____

State: _____ Zip Code: _____ Country: _____

Phone: _____ Attention: _____

Shipping Preference:

Ground 2-Day P.M. 2-Day A.M.

Next Day P.M. Next Day A.M.

(If no preference is indicated, this order will be shipped 2 Day P.M.)


Note: We do not ship new or repaired braces directly to patients.

ROM Post-Op Knee Braces

ROM-SS -- Quantity: _____

Breeze Pads or Wrap-Around Pads

17-22 Inch Length or 19-24 Inch Length



The hinges on these six ROM post-op knee braces feature adjustable flexion and extension stops and a drop lock mechanism that fixates the brace at full extension (an optional -10° locking hinge is available).

ROM-4 -- Quantity: _____

17" Bars (no length adjustment)

Full Shell Pads

- Standard 0° Hinge; or
- Optional -10° Hinge
- Full Shell Black Pads; or
- Full Shell Gray Pads
- Standard Velcro Closures; or
- Optional Buckle Closures

ROM6 -- Quantity: _____

23" Bars (can be shortened to 17")

Full Shell Pads

- Standard 0° Hinge; or
- Optional -10° Hinge
- Full Shell Black Pads; or
- Full Shell Gray Pads
- Standard Velcro Closures; or
- Optional Buckle Closures

ROM6-L -- Quantity: _____

25" Bars (can be shortened to 17")

Full Shell Pads

- Standard 0° Hinge; or
- Optional -10° Hinge
- Full Shell Black Pads; or
- Full Shell Gray Pads
- Standard Velcro Closures; or
- Optional Buckle Closures

AIR ROM-LONG -- Quantity: _____

23" Bars (can be shortened to 17")

AIR Pads

- Standard 0° Hinge; or
- Optional -10° Hinge
- AIR ROM Black Pads; or
- AIR ROM Gray Pads
- Standard Velcro Closures; or
- Optional Buckle Closures

AIR ROM-SHORT -- Quantity: _____

17" Bars (no length adjustment)

AIR Pads

- Standard 0° Hinge; or
- Optional -10° Hinge
- AIR ROM Black Pads; or
- AIR ROM Gray Pads
- Standard Velcro Closures; or
- Optional Buckle Closures

ROM-ES -- Quantity: _____

Full Pads EZ Pads

- Straps With BUCKLE Closures (Full Pads, only) or
- Straps With VELCRO TAB Closures
- STANDARD Length (15-25 Inches) or
- MAXIMUM Length (15-26 Inches)

Special Instructions: _____

Ordering Instructions

Please complete and fax this form to 800.798.2722 (24-hours a day).
 If you are calling in your order, this form indicates the options and information that will be required by our staff. For phone orders, please call 800.700.2722 between 7:00 a.m. and 4:00 p.m. (PST).

Townsend Design

4615 Shepard St., Bakersfield, CA, 93313
 Phone: 800.432.3466 or 661.837.1795; Fax: 800.798.2722