

# Rebel Series

## Ligament Bracing Solutions

### Routinely Prescribed For...

Patients who have experienced ligament injuries and/or reconstructive knee surgery.

### This Product Can Be Ordered...

- Standard OTS sizes XS to XL (XXL Pro model)
- "Customized" from 3 leg measurements
- Custom from a cast, digital scan, tracing, CAS or KneePix

### Offered With These Standard Features...

- Rigid Aircraft Aluminum Frame
- Patented Townsend Motion TM5+ Hinges
- Extension Stop Kit
- Patented Synergistic Suspension Strap
- Anti-Migration Strap Padding
- Rotation Control Tibial Shell Contour
- Powdercoat Finish

### Available Options Include...

- Rebel Lite (lowest weight, less rigid)
- Rebel (ideal for most patients)
- Rebel Pro (heavy duty users and contact sports)
- ACL Or Combined Instabilities
- 13, 14, 15 And 16 Inch Shell Length
- Compression/Suspension Package
- Extension Assist (Pro Only)
- Flexion Stop Kit
- A Variety Of Colors



Townsend's Rebel Series offers an extremely diverse range of ordering options. Based on the patient's clinical needs and activities, you can configure the brace to optimize control and patient compliance.



Optional MC thigh, (available on RebelPro only) provides increased surface contact for added stabilization.

## Exclusive Customized Fabrication

Townsend is the only orthopedic manufacturer that offers the fit quality and value of a custom brace for a pre-sized price when you order any Rebel Series brace from three leg measurements (see page 24 for instructions). For disproportional legs that need to be fit with a non-custom brace, Townsend's "customized" fabrication is the ideal solution. Our staff will mix and match part sizes and pre-set the width between the hinges to minimize or eliminate the need for you to make fitting adjustments.



 **Townsend**

[www.townsenddesign.com](http://www.townsenddesign.com)

Toll Free: 800.432.3466 Ph: 661.837.1795 Fax: 661.837.0613

# General Information

Ordered by : \_\_\_\_\_

Phone: \_\_\_\_\_

Patient's Name: \_\_\_\_\_

BILLING: P.O. Number \_\_\_\_\_

Townsend Account Number: \_\_\_\_\_

Bill To: \_\_\_\_\_

Address: \_\_\_\_\_

City: \_\_\_\_\_

State: \_\_\_\_\_ Zip Code: \_\_\_\_\_ Country: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Ship To: \_\_\_\_\_

Address: \_\_\_\_\_

City: \_\_\_\_\_

State: \_\_\_\_\_ Zip Code: \_\_\_\_\_ Country: \_\_\_\_\_

Phone: \_\_\_\_\_ Attention: \_\_\_\_\_

**Fit Date: If known, please indicate the date you are scheduled to fit the patient:** \_\_\_\_\_

### Shipping Preference:

Ground  2-Day P.M.  2-Day A.M.

Next Day P.M.  Next Day A.M.

(If no preference is indicated, this order will be shipped 2 Day P.M.)

**Note:** We do not ship new or repaired products directly to patients.

# Sizing Reference

### X-Small

6" above: 12.5" to 15.5" 6" below: 11" to 12.5" M-L width: 3" to 3.5"

### Small

6" above: 15.5" to 18.5" 6" below: 12.25" to 13.75" M-L width: 3.5" to 4"

### Medium

6" above: 18.5" to 21" 6" below: 13.25" to 15" M-L width: 4" to 4.5"

### Large

6" above: 21" to 23.5" 6" below: 14.25" to 15.75" M-L width: 4.5" to 5"

### X-Large

6" above: 23.5" to 25" 6" below: 15" to 17" M-L width: 5" to 5.5"

### XX-Large (Pro model only)

6" above: 25" to 28" 6" below: 17" to 19" M-L width: 5.5" to 6"

# Ordering Instructions

Please complete and fax this form to 800.798.2722 (24-hours a day). If you are calling in your order, this form indicates the options and information that will be required by our staff. For phone orders, please call 800.700.2722 between 7:00 a.m. and 4:00 p.m. (PST).

# Rebel Series Ligament Knee Braces

\*Indicates additional charges apply

### Model

#### Rebel

ACL

Combined

#### Rebel Pro

ACL

Combined

#### Rebel Lite

(only available in 13" length)

ACL

Combined

Leg  Right  Left


Thigh Shell Length  7 inch  8 inch

Tibia Shell Length  6 inch  7 inch  8 inch

### Thigh Shell

Standard Thigh (A)

MC Thigh (B) (Pro Model Only)



**TM5+ Hinges** (includes extension stop kit)

Optional Flexion Stop Kit\*

Add optional extension assist bands/posts\* (Pro Model Only)

### Color

Standard colors in stock

Argento Gray  Textured Black

Non-stock colors, brace may ship next business day

Bengal Silver  Bengal Yellow  Bengal White

Gloss Black  Electric Blue  Sky Blue  Violet

Candy Green  Sparkle Copper  Sparkle Red

### CS Package\*

Semi-rigid padded inserts attached inside the thigh shell for dynamic compression and enhanced suspension



### Spooner Patella Stabilizing Attachment\*

### Brace Cover\*

Posterior Closure Style  Pull On Style

Color:  Black  Blue

### Undersleeves\*

18" Cotton  18" Neoprene  22" Neoprene

### Thigh Sleeves\*

1/8 Atrophy Thigh Sleeve  1/16 Comfort Thigh Sleeve

Special Instructions: \_\_\_\_\_


\_\_\_\_\_

\_\_\_\_\_

If your patient has proportional leg sizing (see sizing reference, above) or if you are ordering a brace for stock inventory please select from the size options. However, if you would prefer to have Townsend customize the assembly of your patient's brace at no addition charge please provide leg measurements beside the illustration, below.

Size  X-Small  Small  Medium  Large  
 X-Large  XX-Large (Pro model only)

3 MEASUREMENTS: "CUSTOMIZED"  
ASSEMBLY (NO ADDED CHARGE)



## Townsend Design

4615 Shepard St., Bakersfield, CA, 93313

Phone: 800.432.3466 or 661.837.1795; Fax: 800.798.2722