




The RebelLock is a lightweight option for those patients who need stability and functionality to control significant instabilities.



### Compression Suspension Package

C/S Package enhances brace suspension and increases dynamic support. Semi-rigid padded inserts attached inside the thigh shell compress into the sides of the leg proximal to the knee joints. Anterior and posterior straps allow the patient to adjust the pressure. The C/S Package is a standard feature on all RebelLock knee braces

[www.townsenddesign.com](http://www.townsenddesign.com)

Toll Free: 800.432.3466 Ph: 661.837.1795 Fax: 661.837.0613

# RebelLock

## OTS Specialty Bracing Solutions

### Routinely Prescribed For...

Patients with significant instabilities; weak extensor muscles, failed TKR, limiting or preventing knee flexion.

### This Product Can Be Ordered...

- Standard OTS sizes XS to XL
- "Customized" from 3 leg measurements

### Offered With These Standard Features...

- Aircraft Aluminum Frame
- Drop Lock With Free Motion
- Compression Suspension Package
- Adjustable Locking Position
- Patented Synergistic Suspension Strap
- Anti-Migration Strap Padding
- Rotation Control Tibial Shell Contouring
- Matte And High Gloss Powdercoat Finish

### Available Options Include...

- Shell Length
- Extension Assist Band & Posts

 **Townsend**

## General Information

Ordered by : \_\_\_\_\_

Phone: \_\_\_\_\_

Patient's Name: \_\_\_\_\_

BILLING: P.O. Number \_\_\_\_\_

Townsend Account Number: \_\_\_\_\_

Bill To: \_\_\_\_\_

Address: \_\_\_\_\_

City: \_\_\_\_\_

State: \_\_\_\_\_ Zip Code: \_\_\_\_\_ Country: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Ship To: \_\_\_\_\_

Address: \_\_\_\_\_

City: \_\_\_\_\_

State: \_\_\_\_\_ Zip Code: \_\_\_\_\_ Country: \_\_\_\_\_

Phone: \_\_\_\_\_ Attention: \_\_\_\_\_

**Fit Date: If known, please indicate the date you are scheduled to fit the patient:** \_\_\_\_\_

### Shipping Preference:

Ground  2-Day P.M.  2-Day A.M.

Next Day P.M.  Next Day A.M.

(If no preference is indicated, this order will be shipped 2 Day P.M.)

**Note:** We do not ship new or repaired products directly to patients.

## Sizing Reference

### X-Small

6" above: 12.5" to 15.5" 6" below: 11" to 12.5" M-L width: 3" to 3.5"

### Small

6" above: 15.5" to 18.5" 6" below: 12.25" to 13.75" M-L width: 3.5" to 4"

### Medium

6" above: 18.5" to 21" 6" below: 13.25" to 15" M-L width: 4" to 4.5"

### Large

6" above: 21" to 23.5" 6" below: 14.25" to 15.75" M-L width: 4.5" to 5"

### X-Large

6" above: 23.5" to 25" 6" below: 15" to 17" M-L width: 5" to 5.5"

## Ordering Instructions

Please complete and fax this form to 800.798.2722 (24-hours a day). If you are calling in your order, this form indicates the options and information that will be required by our staff. For phone orders, please call 800.700.2722 between 7:00 a.m. and 4:00 p.m. (PST).

## RebelLock OTS Specialty Knee Braces

\*Indicates additional charges apply

Leg  Right  Left

Thigh Shell Length  7 inch  8 inch

Tibia Shell Length  7 inch  8 inch

Extension Stop Position (factory applied setting)

Zero Degree Terminal Extension/Locking Position

Five Degree Terminal Extension/Locking Position

Hinge Options (includes extension stop kit)

Flexion Stop Kit\*

Add Extension Assist Bands/Posts\*

### Color

Matte Black (standard stock color)

Non-stock colors, brace may ship next business day

Bengal Silver  Bengal Yellow  Bengal White

Gloss Black  Electric Blue  Sky Blue  Violet

Candy Green  Sparkle Copper  Sparkle Red

### Undersleeves\*

18" Cotton  18" Neoprene  22" Neoprene

### Thigh Sleeves\*

1/8 Atrophy Thigh Sleeve  1/16 Comfort Thigh Sleeve

Special Instructions: \_\_\_\_\_


\_\_\_\_\_

\_\_\_\_\_

Size  X-Small  Small  Medium

Large  X-Large

3 MEASUREMENTS: "CUSTOMIZED"  
ASSEMBLY (NO ADDED CHARGE)



## Townsend Design

4615 Shepard St., Bakersfield, CA, 93313

Phone: 800.432.3466 or 661.837.1795; Fax: 800.798.2722