

Rebel Reliever

OA Bracing Solutions

Routinely Prescribed For...

Patients diagnosed with mild to severe medial or lateral OA, as well as patients who have OA in combination with ligament instability.

This Product Can Be Ordered...

- Standard OTS sizes S to XXL (tibia shell available in XS)
- “Customized” from 3 leg measurements
- Custom from a cast, digital scan, tracing, CAS or KneePix

Offered With These Standard Features...

- 60/61 Aircraft Aluminum Frame
- Patented Townsend Motion TM5+ Hinges
- Extension Stop Kit
- Patented Loadshifter
- Patented Synergistic Suspension Strap
- Anti-Migration Strap Padding
- Rotation Control Tibial Shell Contour
- Powdercoat Finish

Available Options Include...


- Medial, Lateral Or Universal
- Shell Length
- Compression/Suspension Package
- PCL Strap
- Flexion Stop Kit
- A Variety Of Colors



The Rebel Reliever is a rigid aircraft aluminum brace designed to help maintain normal leg alignment and shift load away from the damaged side of the knee. These pictures show the optional C/S Package that enhances compression and suspension.

Patented Loadshifter

The brace is fabricated with dual Loadshifter mechanisms for maximum adjustability. This brace can be used universally for medial or lateral OA. Medical professionals can shorten and/or lengthen the uprights to shift the thigh shell angle up to 18 degrees. Increasing corrective force at the proximal end amplifies the affect of the three-point pressure and makes the brace lower profile (at knee center) than other manufacturers' OA braces.



TOWNSEND DESIGN
LOADSHIFTER
RELIEF MECHANISM

 **Townsend**

www.townsenddesign.com

Toll Free: 800.432.3466 Ph: 661.837.1795 Fax: 661.837.0613

General Information

Ordered by : _____

Phone: _____

Patient's Name: _____

BILLING: P.O. Number _____

Townsend Account Number: _____

Bill To: _____

Address: _____

City: _____

State: _____ Zip Code: _____ Country: _____

Phone: _____ Fax: _____

Ship To: _____

Address: _____

City: _____

State: _____ Zip Code: _____ Country: _____

Phone: _____ Attention: _____

Fit Date: If known, please indicate the date you are scheduled to fit the patient: _____

Shipping Preference:

Ground 2-Day P.M. 2-Day A.M.

Next Day P.M. Next Day A.M.

(If no preference is indicated, this order will be shipped 2 Day P.M.)

Note: We do not ship new or repaired products directly to patients.

Sizing Reference

X-Small (tibia shell only)

6" above: 15.5" to 18.5" 6" below: 11" to 12.5" M-L width: 3" to 3.5"

Small

6" above: 15.5" to 18.5" 6" below: 12.25" to 13.75" M-L width: 3.5" to 4"

Medium

6" above: 18.5" to 21" 6" below: 13.25" to 15" M-L width: 4" to 4.5"

Large

6" above: 21" to 23.5" 6" below: 14.25" to 15.75" M-L width: 4.5" to 5"

X-Large

6" above: 23.5" to 25" 6" below: 15" to 17" M-L width: 5" to 5.5"

XX-Large

6" above: 25" to 28" 6" below: 17" to 19" M-L width: 5.5" to 6"

Ordering Instructions

Please complete and fax this form to 800.798.2722 (24-hours a day). If you are calling in your order, this form indicates the options and information that will be required by our staff. For phone orders, please call 800.700.2722 between 7:00 a.m. and 4:00 p.m. (PST).

Rebel Reliever OA Unloading Knee Braces

*Indicates additional charges apply

Compartment

Medial Compartment Lateral Compartment

Universal (select when ordering stock inventory)
Medial OA pads installed, Lateral OA pads included in kit

Leg Right Left

Shell Lengths

Thigh shell extends 7" to 9" on all braces

Tibia 7 inch 8 inch

TM5+ Hinges -- Includes Extension Stops

Optional Flexion Stop Kit*

Color

Standard colors in stock

Argento Gray Textured Black

Non-stock colors, brace may ship next business day


Bengal Silver Bengal Yellow Bengal White

Gloss Black Electric Blue Sky Blue Violet

Candy Green Sparkle Copper Sparkle Red

PCL Strap* (combined instabilities) **Or**

CS Package*
Semi-rigid padded inserts attached inside the thigh shell for dynamic compression and enhanced suspension



Spoooner Patella Stabilizing Attachment*

Brace Cover*

Posterior Closure Style Pull On Style

Color: Black Blue

Undersleeves*

18" Cotton 18" Neoprene 22" Neoprene

Thigh Sleeves*

1/8 Atrophy Thigh Sleeve 1/16 Comfort Thigh Sleeve


Special Instructions: _____

If your patient has proportional leg sizing (see sizing reference, above) or if you are ordering a brace for stock inventory please select from the size options. However, if you would prefer to have Townsend customize the assembly of your patient's brace at no addition charge please provide leg measurements beside the illustration, below.

Size X-Small Small Medium Large

X-Large XX-Large

3 MEASUREMENTS: **"CUSTOMIZED"**
ASSEMBLY (NO ADDED CHARGE)



Townsend Design

4615 Shepard St., Bakersfield, CA, 93313

Phone: 800.432.3466 or 661.837.1795; Fax: 800.798.2722