

SoftForce OA

OA Bracing Solutions

Routinely Prescribed For...

Patients with mild to moderate medial or lateral OA who are more compatible with a less rigid knee brace that offers compression and comfortable leg alignment support.

This Product Can Be Ordered...

- Standard OTS sizes XS to 4XL

Offered With These Standard Features...

- Latex Free Fiber Cool Compression Sleeve
- Patented Townsend Motion TM6 Hinge
- Hinge Pockets (No Exposed Metal)
- Adjustable Aluminum Bar With Paddles
- Condylar Pads
- Extension Stops

Available Options Include...

- Medial Or Lateral OA
- Option Of Lateral Hinge Position For Lateral OA
- Flexion Stop Kit



Townsend's SoftForce provides comfortable correction and inflammation-reducing compression.



The sleeve is designed with a unique dual locking anterior closure thigh that creates a wide opening. This feature makes the brace extremely easy to pull on and tighten.

Hinge/Bar Assembly

The SoftForce can be configured for medial or lateral OA. The aluminum bar and hinge assembly inserts into an internal pocket. For stock inventory, the brace can be used universally by shifting the hinge to the opposite side of the brace.



 **Townsend**

www.townsenddesign.com

Toll Free: 800.432.3466 Ph: 661.837.1795 Fax: 661.837.0613

General Information

Ordered by : _____

Phone: _____

Patient's Name: _____

BILLING: P.O. Number _____

Townsend Account Number: _____

Bill To: _____

Address: _____

City: _____

State: _____ Zip Code: _____ Country: _____

Phone: _____ Fax: _____

Ship To: _____

Address: _____

City: _____

State: _____ Zip Code: _____ Country: _____

Phone: _____ Attention: _____

Fit Date: If known, please indicate the date you are scheduled to fit the patient: _____

Shipping Preference:

Ground 2-Day P.M. 2-Day A.M.

Next Day P.M. Next Day A.M.

(If no preference is indicated, this brace will be shipped 2 Day P.M.)

Note: We do not ship new or repaired braces directly to patients.

Sizing Reference

With the patient's leg at full extension, measure circumference around the center of the affected knee. Also, reference minimum and maximum calf and thigh circumferences to determine the best size.

X-Small

12" to 13" knee circumference (calf 12" to 13" and thigh 16" to 24")

Small

13" to 14" knee circumference (calf 13" to 15" and thigh 17" to 26")

Medium

14" to 15" knee circumference (calf 14" to 18" and thigh 18" to 28")

Large

15" to 17" knee circumference (calf 15" to 20" and thigh 20" to 29")

X-Large

17" to 19" knee circumference (calf 16" to 22" and thigh 22" to 31")

XX-Large

19" to 21" knee circumference (calf 18" to 24" and thigh 23" to 32")

XXX-Large

21" to 22" knee circumference (calf 20" to 26" and thigh 25" to 34")

XXXX-Large

22" to 23" knee circumference (calf 21" to 27" and thigh 26" to 36")

SoftForce OA Unloading Knee Braces

***Indicates additional charges apply**

Compartment Medial OA Lateral OA

If you select lateral OA, please read these instructions and select the lateral OA hinge position.

Optimum unloading for Lateral OA patients is achieved when the hinge is placed on the medial side of the knee. However, if the patient is extremely knock-knee and a medial hinge is problematic, request the Lateral Hinge position.

Medial Hinge Position Lateral Hinge Position

Leg Right Left

Size X-Small Small Medium Large

X-Large XX-Large XXX-Large XXXX-Large

TM6 Hinge -- Includes extension stop kit

Optional Flexion Stop Kit*

Correction

(hinge bar will be pre-bent)

3° Standard Corrective Force

5° Moderate Corrective Force

7° Maximum Corrective Force

Replacement Sleeves

Indicate Replacement Sleeve Size

X-Small Small Medium

Large X-Large XX-Large

XXX-Large XXXX-Large

Ordering Instructions

Please complete and fax this form to 800.798.2722 (24-hours a day). If you are calling in your order, this form indicates the options and information that will be required by our staff. For phone orders, please call 800.700.2722 between 7:00 a.m. and 4:00 p.m. (PST).

For Multiple Brace (Stock Inventory) Orders

| Quantity | Brace Size | For stock braces, unless otherwise specified, the bars are pre-bent for 3 degrees of correction. Sleeves are universal (left or right leg). As such, a sleeve with a bar inserted in the pocket on the right side of the sleeve can be used as a right leg medial OA brace, or as a left leg lateral OA brace (with a medial hinge). A sleeve with a bar inserted in the pocket on the left side of the sleeve can be used as a left leg medial OA brace, or as a right leg lateral OA brace (with a medial hinge). If you need to use a brace for lateral OA, with a lateral hinge position, the bar must be re-bent. Instructions are provided with the brace to ensure the fitter has the hinge/ bar assembly inserted into the appropriate side of the sleeve and bent in a direction that will provide the desired unloading forces. When ordering, you can request that we fill orders for a specific leg and compartment. But to minimize inventory, you can maintain one brace in each size and reposition (and potentially re-bend) the bar as needed. |
|--|------------|---|
| | X-Small | |
| | Small | |
| | Medium | |
| | Large | |
| | X-Large | |
| | XX-Large | |
| | XXX-Large | |
| | XXXX-Large | |
| NOTE: Braces shipped for stock inventory will be packaged with the hinge inserted in the right pocket (right leg medial, left leg lateral). | | |

Townsend Design

4615 Shepard St., Bakersfield, CA, 93313

Phone: 800.432.3466 or 661.837.1795; Fax: 800.798.2722