

Distal Extension To Townsend Custom Knee Braces

Complete the information below and submit this form with the **CUSTOM BRACE ORDER FORM** you fill out for this patient.

Patient Last Name _____

Patient First Name _____

Brace Model _____

Account Name _____

REQUIRED MEASUREMENTS


Please Take Measurements For Placement Of Knee Hinges & Ankle Joint

Measurement from floor

to knee center is _____ inches.

Measurement from floor to the lateral

apex of the ankle is _____ inches.



SUSPENSION/ROTATION EXTENSION

(Note: 2-3 additional days of production are required for braces with a distal extension)

Extension Options

- Extend on Lateral Side of Brace (standard)
 Extend on Medial Side of Brace

Foot Plate Options

Choose Material:


- Polypro (adjustable) Graphite (nonadjustable, thicker)

Choose Length:

- Heel Cup Trim to Toe Sulcus Trim Behind Met. Heads

Choose Heel Option:

- Open Heel Cup (graphite, only) 1/2 Heel Full Heel



The extension can be added to Premier Series (left) and traditional full shell braces (right). The option is available on functional or OA braces designs.

Description

The extension includes a graphite bar laminated to the tibia shell; a free ankle joint; and a graphite or polypro foot plate. The primary function of the extension is to prevent migration and/or provide additional rotation control.

Casting

A full leg cast is required including the foot (weight-bearing).

Ordering

In addition to filling out this form, you must also submit the appropriate order form for the type of brace you are ordering (Premier, Air Townsend, etc.) as this form does not indicate patient information, billing and shipping information, brace model, color and options on the order form. An extension cannot be added to an existing brace. A new brace must be fabricated with the extension.

Updated: June, 2006