

General Information

Ordered by : _____

Phone: _____

BILLING: P.O. Number _____

Townsend Account Number: _____

Bill To: _____

Address: _____

City: _____

State: _____ Zip Code: _____ Country: _____

Phone: _____ Fax: _____

Ship To: _____

Address: _____

City: _____

State: _____ Zip Code: _____ Country: _____

Phone: _____ Attention: _____

Fit Date: If known, please indicate the date you are scheduled to fit the patient: _____

Shipping Preference:

Ground 2-Day P.M. 2-Day A.M.

Next Day P.M. Next Day A.M.

(If no preference is indicated, this order will be shipped 2 Day P.M.)

Note: We do not ship new or repaired products directly to patients.

Special Instructions: _____

FAB Series & FG Series Casting Tape

FAB SERIES

This is a unique elastic polyester blend casting tape (not fiber glass, not plaster) that is easier to contour than traditional fiber glass. It has the molding characteristics of plaster without the mess.

FAB2 (2" X 4 yards; white tape)

Boxes: 1 2 3 4 Other _____

FAB3 (3" X 4 yards; white tape)

Boxes: 1 2 3 4 Other _____

FAB4 (4" X 4 yards; white tape)

Boxes: 1 2 3 4 Other _____

FAB5 (5" X 4 yards; white tape)

Boxes: 1 2 3 4 Other _____

FG SERIES

This is a multi-purpose fiber glass casting tape for use/wear by a patient for stabilization or positional support following an injury or surgery. The FG series is also suitable for creating patient molds for use in the fabrication of orthopedic braces and prosthetic devices.

FG2 (2" X 4 yards; cast tape)

Color: White Red Lime Green

Boxes: 1 2 3 4 Other _____

FG3 (3" X 4 yards; cast tape)

Color: White Red Lime Green

Boxes: 1 2 3 4 Other _____

FG4 (4" X 4 yards; cast tape)

Color: White Red Lime Green

Boxes: 1 2 3 4 Other _____


FG5 (5" X 4 yards; white tape)

Boxes: 1 2 3 4 Other _____


Instructions For Use

NOTE: KEEP TAPE WET


Tear Top Off Of Bag, Leave Roll In Bag, & Immerse In Warm Water



Squeeze Out Excess Water (With Roll In Bag), Submerge/ Squeeze Again



Remove Roll From Bag, Roll Onto Limb & Contour



Allow To Dry Until Rigid & Use Shark-Nose Scissors To Remove

TOWNSEND'S SHIPPING DEPARTMENT USE ONLY

NEW BRACE SERVICE -- Original Brace Returned? Yes No

RECEIVED _____ SHIPPED VIA _____

Townsend Design

4615 Shepard St., Bakersfield, CA, 93313
Phone: 800.432.3466 or 661.837.1795; Fax: 800.798.2722