

General Information

Ordered by : _____

Phone: _____

Patient's Name: _____

BILLING: P.O. Number _____

Townsend Account Number: _____

Bill To: _____

Address: _____

City: _____

State: _____ Zip Code: _____ Country: _____

Phone: _____ Fax: _____

Ship To: _____

Address: _____

City: _____

State: _____ Zip Code: _____ Country: _____

Phone: _____ Attention: _____

Fit Date: If known, please indicate the date you are scheduled to fit the patient: _____

Shipping Preference:

Ground 2-Day P.M. 2-Day A.M.

Next Day P.M. Next Day A.M.

(If no preference is indicated, this order will be shipped 2 Day P.M.)

Note: We do not ship new or repaired products directly to patients.

Sizing Reference

Measure from the wrist up to the elbow crease. If this measurement is eight inches or less, then a Small/Medium size or 15 inch model should be ordered. If the measurement is longer than eight inches, then it is best to order the Large/X-Large size or 16 inch model.

Ordering Instructions

Please complete and fax this form to 800.798.2722 (24-hours a day). If you are calling in your order, this form indicates the options and information that will be required by our staff. For phone orders, please call 800.700.2722 between 7:00 a.m. and 4:00 p.m. (PST).

ROM Post-Op Elbow Braces

ROM-LBO

The LBO comes standard with four straps, tool less adjustability, 8 individual Breeze pads or 4 individual full wrap pads; single or dual upright applications. Easy adjustable flexion and extension stops provide range of motion control in 10 degree increments. The LBO can be locked at 90 degrees for post-surgical protocols. For single hinge applications, right or left arm must be indicated.

Single Hinge (SH) Model

Quantity: ___ LEFT ARM

Quantity: ___ RIGHT ARM

Pad Options:

Breeze Pads

Individual Full Wrap

Dual Hinge (DH) Model

Quantity: ___ UNIVERSAL

Pad Options:

Breeze Pads

Individual Full Wrap



Multi-Adjustment ROM Elbow

This brace has four straps and individual black pads that fit under each strap. Single hinge/bar and dual hinge/bar models can be ordered with a 15 or 16 inch bar length. Flexion and extension stops are in 10 degree increments and can be set without tools. Locking is possible by offsetting the flexion and extension stops at the same degree. For single hinge braces, you must indicate right or left arm.

Single Hinge (SH) Model

Quantity: ___ LEFT ARM Length: ___ 15 inch ___ 16 inch

Quantity: ___ RIGHT ARM Length: ___ 15 inch ___ 16 inch

Dual Hinge (DH) Model

Length (Indicate quantity): ___ 15 inch ___ 16 inch

Townsend's Multi-Adjustment ROM Elbow Brace features adjustable flexion and extension stops.



BRACES WITH LOW PROFILE HINGES

Flexion can be limited at 90, 110 and 120 degrees, and extension can be stopped at 0, 15, 30, 45 and 90 degrees. The screw stops are adjusted using a small flathead screwdriver. The hinge can be locked at 90 degrees.


Low Profile ROM Elbow (Dual Hinges) With Shells And Optional Hand Extension

The DUAL HINGE (DH) model is fabricated with rigid ABS plastic shells. This brace can be ordered with an optional bar and hand extension for increased stabilization and rotation control. The shells conform to the arm when the straps are tightened. Also, the length of the hand extension and position of the shells is adjustable.

Quantity: ___ S/M ___ L/XL

Install Optional Hand Extension? ___ Yes ___ No

Include Shoulder Sling Strap?* ___ Yes ___ No



Townsend Design

4615 Shepard St., Bakersfield, CA, 93313
Phone: 800.432.3466 or 661.837.1795; Fax: 800.798.2722

Special Instructions:
