

General Information

Ordered by : _____

Phone: _____

Patient's Name: _____

BILLING: P.O. Number _____

Townsend Account Number: _____

Bill To: _____

Address: _____

City: _____

State: _____ Zip Code: _____ Country: _____

Phone: _____ Fax: _____

Ship To: _____

Address: _____

City: _____

State: _____ Zip Code: _____ Country: _____

Phone: _____ Attention: _____

Fit Date: If known, please indicate the date you are scheduled to fit the patient: _____

Shipping Preference:

Ground 2-Day P.M. 2-Day A.M.

Next Day P.M. Next Day A.M.

(If no preference is indicated, this order will be shipped 2 Day P.M.)

Note: We do not ship new or repaired products directly to patients.

Sizing Reference

X-Small

6" above: 12.5" to 15.5" 6" below: 11" to 12.5" M-L width: 3" to 3.5"

Small

6" above: 15.5" to 18.5" 6" below: 12.25" to 13.75" M-L width: 3.5" to 4"

Medium

6" above: 18.5" to 21" 6" below: 13.25" to 15" M-L width: 4" to 4.5"

Large

6" above: 21" to 23.5" 6" below: 14.25" to 15.75" M-L width: 4.5" to 5"

X-Large

6" above: 23.5" to 25" 6" below: 15" to 17" M-L width: 5" to 5.5"

XX-Large (Pro model only)

6" above: 25" to 28" 6" below: 17" to 19" M-L width: 5.5" to 6"

Ordering Instructions

Please complete and fax this form to 800.798.2722 (24-hours a day). If you are calling in your order, this form indicates the options and information that will be required by our staff. For phone orders, please call 800.700.2722 between 7:00 a.m. and 4:00 p.m. (PST).

Rebel Series Ligament Knee Braces

*Indicates additional charges apply

Model

Rebel

ACL

Combined

Rebel Pro

ACL

Combined

Rebel Lite

(only available in 13" length)

ACL

Combined

Leg Right Left


Thigh Shell Length 7 inch 8 inch

Tibia Shell Length 6 inch 7 inch 8 inch

Thigh Shell

Standard Thigh (A)

MC Thigh (B) (Pro Model Only)



TM5+ Hinges (includes extension stop kit)

Optional Flexion Stop Kit*

Add optional extension assist bands/posts* (Pro Model Only)

Color

Standard colors in stock

Argento Gray Textured Black

Non-stock colors, brace may ship next business day

Bengal Silver Bengal Yellow Bengal White

Gloss Black Electric Blue Sky Blue Violet

Candy Green Sparkle Copper Sparkle Red

CS Package*

Semi-rigid padded inserts attached inside the thigh shell for dynamic compression and enhanced suspension



Spooner Patella Stabilizing Attachment*

Brace Cover*

Posterior Closure Style Pull On Style

Color: Black Blue

Undersleeves*

18" Cotton 18" Neoprene 22" Neoprene

Thigh Sleeves*


1/8 Atrophy Thigh Sleeve 1/16 Comfort Thigh Sleeve

Special Instructions: _____

If your patient has proportional leg sizing (see sizing reference, above) or if you are ordering a brace for stock inventory please select from the size options. However, if you would prefer to have Townsend customize the assembly of your patient's brace at no addition charge please provide leg measurements beside the illustration, below.

Size X-Small Small Medium Large
 X-Large XX-Large (Pro model only)

3 MEASUREMENTS: "CUSTOMIZED"
ASSEMBLY (NO ADDED CHARGE)



Townsend Design

4615 Shepard St., Bakersfield, CA, 93313

Phone: 800.432.3466 or 661.837.1795; Fax: 800.798.2722